

1701 AR•Model PCH Articulator

1610 AR•Magnetic Model PCH Articulator

The model PCH articulator is a precision articulator which allows the transfer of dental casts from one PCH or PSH articulator to another. The model PCH has a curved pin calibrated in degrees.

The 1610 AR Comes with 20 Magnetic Mounting Plates & Hex Wrench.



1801 AR•Model PSH Articulator

1620 AR•Magnetic Model PSH Articulator

The model PSH articulator is a precision articulator which allows the transfer of dental casts from one PSH or PCH articulator to another. The model PSH comes with a straight incisal pin calibrated in millimeters.

The 1620 AR Comes with 20 Magnetic Mounting Plates & Hex Wrench.



1650 AR•Magnetic Model PCF Articulator

The model PCF articulator is a precision articulator featuring fixed 1.5mm analogs set at 20°, an installed magnetic system and a curved incisal pin calibrated in degrees.

The 1650 AR Comes with 10 Magnetic Mounting Plates.

