



4000 FB•Pana-Mount™ Face-Bow

The Pana-Mount™ Face-Bow features wrench or wrenchless detachable bite-fork assemblies which allows the face-bow to be used for multiple patients without having to mount casts immediately. It has a standard 22mm nasion relator as well as an orbital pointer. Comfortable ear pieces and foam nasion pad increase patient comfort.

Includes:

Face-Bow, two complete Bite-Fork Assemblies, Mounting Fixture, Standard Nasion Relator, T-Screws, and Hex Wrench



4070 FB•Bio-Esthetic™ Level Gauge

The Bio-Esthetic™ Level Gauge is used to level the face-bow horizontal so the anterior esthetic plane of the teeth can be established correctly on the articulator.



4027 FB•Bite-Fork Assembly (Complete)

The Bite-Fork Assembly is used with the Pana-Mount™ Face-Bow and comes complete with a standard bite-fork and stem assembly, as well as T-Screws.



4028 FB•Bite-Fork (Regular)

The standard bite-fork is used with the Pana-Mount™ Face-Bow and bite-fork stem assembly. Constructed of all stainless steel, it is designed for easy cleaning and sterilization.

4037 FB•Bite-Fork (Full denture)

The full denture bite-fork is used with the Pana-Mount™ Face-Bow and bite-fork stem assembly for transferring maxillary bite-rim and edentulous cast to the articulator. Constructed of all stainless steel, it is designed for easy cleaning and sterilization.

4038 FB•Bite-Fork (Pedo)

The pedo bite-fork is used with the Pana-Mount™ Face-Bow and bite-fork stem assembly. Constructed of all stainless steel, it is designed for easy cleaning and sterilization.



4080 FB•Adjustable Nasion Relator

The Adjustable Nasion Relator is used in conjunction with the Bio-Esthetic™ Level Gauge to level the face-bow in the sagittal plane when the patient's head is erect. The adjustable nasion is calibrated in millimeters.



4030 FB•Bite-Fork Stem Assembly (Only)

This stem assembly is for attaching the bite-fork to the Pana-Mount™ Face-Bow. It is the same assembly as item #4027, excluding the bite-fork. The T-Screws are also included.