

4300 KS•Kois Dento-Facial Analyzer System

The Kois Dento-Facial Analyzer System registers and transfers the patient's occlusal plane in three planes of space to the articulator related to an average 100mm axis-incisal distance.

Includes:

#4305 KS•Kois Dento-Facial Analyzer

#4360 KS•Kois Adjustable Platform




4305 FB•Kois Dento-Facial Analyzer

The Kois Dento-Facial Analyzer registers and transfers the patient's occlusal plane in three planes of space to the articulator.

Includes:

#4322 KS•Kois disposable Index Trays (20)



4322 KS•Kois Disposable Index Trays (20)

4323 KS•Kois Disposable Index Trays (50)

4324 KS•Kois Disposable Index Trays (100)

The disposable index trays are used with the Kois Dento-Facial Analyzer and provides a permanent or disposable bite-fork registration record.



4360 KS•Kois Adjustable Platform

The Kois Adjustable Platform is used with the Kois Index Trays to mount the patient's maxillary cast to the articulator. The platform can be adjusted up or down to a new incisal edge or occlusal plane position.

Includes:

4331 Standard Waxing Guide

