



2850 MP • Magna-Split™ II System

This magnetic mounting plate system is used with the PCH or PSH articulator. It comes with removable rare earth magnet and stainless steel screws for fixing the plate to the upper and lower frames of the articulator.

Includes:

2-2804 MP • Magna-Split™ Mounting Plate

1-2860 MP • Magna-Split™ Index Plates (20)



2860 MP • Magna-Split™ Index Plates (20)

2861 MP • Magna-Split™ Index Plates (50)

2862 MP • Magna-Split™ Index Plates (100)

These pre-formed plastic Index Plates attach directly to the Magna-Split™ Mounting Plate for easy mounting of study casts and can be used with our PAL™ - Plastic Articulator Lite .



2804 MP • Magna-Split™ Mounting Plate

The mounting plate comes with removable rare earth magnet and two stainless steel screws for fixing the plate to upper or lower frame of the articulator or Axis Mounting System.

(Picture shows magnet in place.)



2852 MP • Magna-Split™ Index Plates-Old Style (20)

2853 MP • Magna-Split™ Index Plates-Old Style (50)

2854 MP • Magna-Split™ Index Plates-Old Style (100)

These pre-formed plastic Index Plates attach directly to the Magna-Split™ Mounting Plate for quick and easy mounting of study casts.



2706 MP • Plastic Mounting Plates (10)

2707 MP • Plastic Mounting Plates (20)

2708 MP • Plastic Mounting Plates (50)

These Plastic Mounting Plates are standard screw on mounting plates. They can be used with several brands of articulators.